

LANE
STAMP

NOV 8 1989

STANFORD, CA 94305

Antibodies


A LABORATORY MANUAL

Ed Harlow

Cold Spring Harbor Laboratory

David Lane

Imperial Cancer Research Fund Laboratories



Cold Spring Harbor Laboratory
1988

LANE BEST AVAILABLE COPY

Antibodies A LABORATORY MANUAL

All rights reserved

© 1988 by Cold Spring Harbor Laboratory

Printed in the United States of America

Book and cover design by Emily Harste

Cover: "Nature Abstracted," watercolor by Carl Molno

Library of Congress Cataloging-in-Publication Data

Antibodies : a laboratory manual / by Ed Harlow, David Lane.

p. cm.

Bibliography: p.

Includes index.

ISBN 0-87969-314-2

1. Immunoglobulins--Laboratory manuals. 2. Immunochemistry--Laboratory manuals. I. Harlow, Ed. II. Lane, David (David P.).

1952--

QR186.7.A53 1988

574.2'93'028--dc19

Researchers using the procedures of this manual do so at their own risk. Cold Spring Harbor Laboratory makes no representations or warranties with respect to the material set forth in this manual and has no liability in connection with the use of these materials.

Certain experimental procedures in this manual may be the subject of national or local legislation or agency restrictions. Users of this manual are responsible for obtaining the relevant permissions, certificates, or licenses in these cases. Neither the authors of this manual nor Cold Spring Harbor Laboratory assume any responsibility for failure of a user to do so.

Authorization to photocopy items for internal or personal use, or the internal or personal use of specific clients, is granted by Cold Spring Harbor Laboratory for libraries and other users registered with the Copyright Clearance Center (CCC) Transactional Reporting Service, provided that the base fee of \$1.00 per article is paid directly to CCC, 27 Congress St., Salem MA 05970. [0-87969-314-2/88 \$1.00 + .00]. This consent does not extend to other kinds of copying, such as copying for general distribution, for advertising or promotional purposes, for creating new collective works, or for resale.

All Cold Spring Harbor Laboratory publications may be ordered directly from Cold Spring Harbor Laboratory, Box 100, Cold Spring Harbor, New York 11724. Phone: 1-800-843-4388. In New York (516) 367-8423.

BEST AVAILABLE COPY

■ PURIFYING ANTIBODIES

Purified antibodies are required for a number of techniques. Table 8.1 lists several techniques that rely on purified antibodies, at least in some steps. In many of the examples listed in this table, purified antibodies are labeled with an easily detected "tag" (p. 319), and these labeled antibodies are then used to determine the presence of an antigen or another antibody. When labeled anti-immunoglobulin antibodies (p. 622) are used to measure the presence of other antibodies, it is seldom worthwhile to prepare and label these reagents yourself. They can be purchased from a number of commercial sources, where they are prepared and tested in large, economic batches. However, when labeled antibodies will be used to detect an antigen directly, the primary antibody must be purified first. Direct labeling also allows two antibodies to be compared in the same assay by marking them with different tags. In other instances, purified antibodies are necessary for different applications. For example, purified antibodies may lower the background in some assays or purification may be the easiest method to concentrate antibody solutions.

There are a wide variety of methods used to purify antibodies. The correct choice of purification method will depend on a number of variables, including the use for which the antibodies are intended, the species in which it was raised, its class and subclass if it is a monoclonal antibody, and the source that will serve as the starting material for the purification. Table 8.2 summarizes the possible sources of antibodies for purification. Also included in this table are the possible sources of antibody contamination and the expected level of purity. Both of these factors may influence the choice of starting material.

TABLE 8.1
Techniques That Require Purified Antibodies

Technique	Antibody use	Antibody type	Best sources	Comments
Cell Staining	Direct localization	Anti-antigen	Polyclonal or monoclonal	Prepare yourself
	Indirect localization	Anti-antibody	Polyclonal	Available commercially
Immunoassays	Direct detection	Anti-antigen	Monoclonal, but polyclonal OK in some cases	Prepare yourself
	Indirect detection	Anti-antibody	Polyclonal	Available commercially
Immunoblots	Direct detection	Anti-antigen	Polyclonal or monoclonal	Prepare yourself
	Indirect detection	Anti-antibody	Polyclonal	Available commercially
Immunoaffinity	Purification	Anti-antigen	Monoclonal	Prepare yourself

BEST AVAILABLE COPY

Table 8.3 summarizes the commonly used methods for antibody purification and also lists the advantages and disadvantages of each method. As can be seen from this table, it is often necessary to combine several methods to achieve the desired purification. Finally, Tables 8.4 and 8.5 compare the expected results of the different purification methods when using either polyclonal (Table 8.4) or monoclonal (Table 8.5) antibody sources.

Although no single suggestion can fulfill all of the requirements for different purification needs, most workers have found that purification on protein A beads is the most useful technique. This technique has become the method of choice for antibodies with high affinities for protein A.

■ Conventional Methods

Antibodies from serum or ascites can be purified using conventional methods involving precipitation and column chromatography. When similar techniques are used on tissue culture supernatants, the degree of purity achieved will be lower because of the lower concentrations of the specific antibodies. For tissue culture supernatants, purification using protein A beads (p. 309) or anti-immunoglobulin antibody affinity columns (p. 316) is recommended.

BEST AVAILABLE COPY